

Charger 24V SLA (MSP-C-AH-00100)

Specifications

Input voltage	100-240 VAC 50-60Hz
Input current	700-500mA
Output voltage	24 Volt DC
Output current	1000mA
IP Class charger	IP40
Certification	UL2601-1 IEC60601-1 IEC60335-2-29
Operational temperature	0°C to +40°C
Storage temperature	-40°C to +70°C
Case colour	Blue RAL5007, Black charger
Dimensions LxWxH	855x110x100mm



Order Information

Article number	MSP-C-AH-00101 for EU MSP-C-AH-00102 for UK MSP-C-AH-00103 for AUS MSP-C-AH-00104 for USA
Compatible with Arjo number	KPA8200 SPL3024
Packaging	Individual packed, 50 pcs per pallet
M.O.Q.	1 piece
Delivery time	2-3 working days

Features

- + Suitable for MSP-A-AH-00200 / Arjo KPA0100
- + Charger specially designed for medical equipment
- + Delivered with connection cable suited for UK, USA, EU or AUS
- + When the battery pack is charged, the charger will continue as a trickle charger

Charger 24V SLA (MSP-C-AH-00200)

Specifications

Input voltage	100-240 VAC 50-60Hz
Input current	700-500mA
Output voltage	24 Volt DC
Output current	1000mA
IP Class charger	IP40
Certification	UL2601-1 IEC60601-1 IEC60335-2-29
Operational temperature	0°C to +40°C
Storage temperature	-40°C to +70°C
Case colour	Blue RAL5007, Black charger
Dimensions LxWxH	715x100x45mm



Order Information

Article number	MSP-C-AH-00201 for EU MSP-C-AH-00202 for UK MSP-C-AH-00203 for AUS MSP-C-AH-00204 for AUS
Compatible with Arjo number	NDA8200 / NDA2200 / HMX987(S)
Packaging	Individual packed, 140 pcs per pallet
M.O.Q.	1 piece
Delivery time	2-3 working days

Features

- + Suitable for MSP-A-AH-001X0 / Arjo NDA0100-XX HMX991 NDA0200
- + Charger specially designed for medical equipment
- + Delivered with connection cable suited for UK, USA, EU or AUS
- + When the battery pack is charged, the charger will continue as a trickle charger

Charger 24V SLA (MSP-C-AH-00300)

Specifications

Input voltage	100-240 VAC 50-60Hz
Input current	700-500mA
Output voltage	24 Volt DC
Output current	1000mA
IP Class charger	IP40
Certification	UL2601-1 IEC60601-1 IEC60335-2-29
Operational temperature	0°C to +40°C
Storage temperature	-40°C to +70°C
Case colour	Blue RAL5007, Black charger
Dimensions LxWxH	405x99x160mm



Order Information

Article number	MSP-C-AH-00301 for EU MSP-C-AH-00302 for UK MSP-C-AH-00303 for AUS MSP-C-AH-00304 for USA
Compatible with Arjo number	NDA1000
Packaging	Individual packed, 50 pcs per pallet
M.O.Q.	1 piece
Delivery time	2-3 working days

Features

- + Suitable for MSP-A-AH-00400 / Arjo KKA1100-XX
- + Charger specially designed for medical equipment
- + Delivered with connection cable suited for UK, USA, EU or AUS
- + When the battery pack is fully charged, the charger will continue as a trickle charger

Charger 24V SLA (MSP-C-AH-00400)

Specifications

Input voltage	100-240 VAC 50-60Hz
Input current	700-500mA
Output voltage	24 Volt DC
Output current	1000mA
IP Class charger	IP40
Certification	UL2601-1 IEC60601-1 IEC60335-2-29
Operational temperature	0°C to +40°C
Storage temperature	-40°C to +70°C
Case colour	Blue RAL5007, Black charger
Dimensions LxWxH	330x95x80mm



Order Information

Article number	MSP-C-AH-00401 for EU MSP-C-AH-00402 for UK MSP-C-AH-00403 for AUS MSP-C-AH-00404 for USA
Compatible with Arjo number	KTA0101
Packaging	Individual packed, 100 pcs per pallet
M.O.Q.	1 piece
Delivery time	2-3 working days

Features

- + Suitable for MSP-A-AH-00700 / Arjo KTA0102
- + Charger specially designed for medical equipment
- + Delivered with connection cable suited for UK, USA, EU or AUS
- + When the battery pack is fully charged, the charger will continue as a trickle charger

Charger 24V/DC NiMH (MSP-C-AH-00500)

Specifications

Input voltage	100-240 VAC 50-60Hz
Input current	700-500mA
Output voltage	24 V/DC
Output current	400mA
IP class charger	IP40
Certification	UL2601-1 IEC60601-1 IEC60335-2-29
Operational temperature	0°C to +40°C
Storage temperature	-40°C to +70°C
Case colour	Black
Dimensions LxWxH	95x48x33mm

Order Information

Article number	MSP-C-AH-00501 for EU MSP-C-AH-00502 for GB MSP-C-AH-00503 for AUS MSP-C-AH-00504 for USA
Compatible with Arjo art. number	NEA1000
Packaging	Individually packed, 50pcs per box
M.O.Q.	1 piece
Delivery Time	2 working days



Features

- + Suitable for MSP-A-AH-00300 / Arjo NEA0100-083 NEA0100
- + Charger specifically designed for medical equipment
- + Delivered with power cord suited for UK, EU, USA or AUS
- + When the battery pack is charged, the charger will continue as trickle charger

Charger 24V/DC SLA (MSP-C-AH-00600)

Specifications

Input voltage	100-240 VAC 50-60Hz
Input current	700-500mA
Output voltage	24 V/DC
Output current	400mA
IP class charger	IP40
Certification	UL2601-1 IEC60601-1 IEC60335-2-29
Operational temperature	0°C to +40°C
Storage temperature	-40°C to +70°C
Case colour	Black
Dimensions LxWxH	109x69x39mm

Order Information

Article number	MSP-C-AH-00601 for EU MSP-C-AH-00602 for GB MSP-C-AH-00603 for AUS MSP-C-AH-00604 for USA
Replacement for Arjo art. number	NDA8200/NDA1200/HMX987
Packaging	Individually packed, 50pcs per box
M.O.Q.	1 piece
Delivery Time	2 working days



Features

- + Suitable for MSP-A-AH-001X0 / Arjo NDA0100/HMX991/NDA0200
- + Charger specifically designed for medical equipment
- + Delivered with power cord suited for UK, EU, USA or AUS
- + When the battery pack is charged, the charger will continue as trickle charger

Charger 24V/DC SLA (MSP-C-AH-00700)

Specifications

Input voltage	100-240 VAC 50-60Hz
Input current	700-500mA
Output voltage	24 V/DC
Output current	400mA
IP class charger	IP40
Certification	UL2601-1 IEC60601-1 IEC60335-2-29
Operational temperature	0°C to +40°C
Storage temperature	-40°C to +70°C
Case colour	Black
Dimensions LxWxH	109x69x39mm

Order Information

Article number	MSP-C-AH-00701 for EU MSP-C-AH-00702 for GB MSP-C-AH-00703 for AUS MSP-C-AH-00704 for USA
Replacement for Arjo art. number	New product
Packaging	Individually packed, 50pcs per box
M.O.Q.	1 piece
Delivery Time	2 working days



Features

- + Suitable for MSP-A-AH-00400 / Arjo KKA1100-04 KKA1100-XX
- + Charger specifically designed for medical equipment
- + Delivered with power cord suited for UK, EU, USA or AUS
- + When the battery pack is charged, the charger will continue as trickle charger

Charger 24V/DC SLA (MSP-C-AH-00800)

Specifications

Input voltage	100-240 VAC 50-60Hz
Input current	700-500mA
Output voltage	24 V/DC
Output current	400mA
IP class charger	IP40
Certification	UL2601-1 IEC60601-1 IEC60335-2-29
Operational temperature	0°C to +40°C
Storage temperature	-40°C to +70°C
Case colour	Black
Dimensions LxWxH	109x69x39mm

Order Information

Article number	MSP-C-AH-00801 for EU MSP-C-AH-00802 for GB MSP-C-AH-00803 for AUS MSP-C-AH-00804 for USA
Replacement for Arjo art. number	KPA2000/KPA8200/SPL3024
Packaging	Individually packed, 50pcs per box
M.O.Q.	1 piece
Delivery Time	2 working days



Features

- + Suitable for MSP-A-AH-00200/500 / Arjo KPA0100/SPL3021/KTA0102
- + Charger specifically designed for medical equipment
- + Delivered with power cord suited for UK, EU, USA or AUS
- + When the battery pack is charged, the charger will continue as trickle charger

Charger 12V/DC SLA (MSP-C-AH-00900)

Specifications

Input voltage	100-240 VAC 50-60Hz
Input current	700-500mA
Output voltage	12 V/DC
Output current	500mA
IP class charger	IP40
Certification	UL2601-1 IEC60601-1 IEC60335-2-29
Operational temperature	0°C to +40°C
Storage temperature	-40°C to +70°C
Case colour	Black
Dimensions LxWxH	405x99x160mm

Order Information

Article number	MSP-C-AH-00901 for EU MSP-C-AH-00902 for GB MSP-C-AH-00903 for AUS MSP-C-AH-00904 for USA
Replacement for Arjo art. number	NCA1000
Packaging	Individually packed, 50pcs per box
M.O.Q.	1 piece
Delivery Time	2 working days



Features

- + Suitable for Concerto / Arjo 102193
- + Charger specifically designed for medical equipment
- + Delivered with power adapter suited for UK, EU, USA or AUS
- + When the battery pack is charged, the charger will continue as trickle charger

Charger 24V/DC SLA (MSP-C-AH-01000)

Specifications

Input voltage	100-240 VAC 50-60Hz
Input current	700-500mA
Output voltage	24 V/DC
Output current	400mA
IP class charger	IP40
Certification	UL2601-1 IEC60601-1 IEC60335-2-29
Operational temperature	0°C to +40°C
Storage temperature	-40°C to +70°C
Case colour	Black
Dimensions LxWxH	109x69x39mm

Order Information

Article number	MSP-C-AH-01001 for EU MSP-C-AH-01002 for GB MSP-C-AH-01003 for AUS MSP-C-AH-01004 for USA
Replacement for Arjo art. number	700-24201
Packaging	Individually packed, 50pcs per box
M.O.Q.	1 piece
Delivery Time	2 working days



Features

- + Suitable for Maxi Sky 440
- + Charger specifically designed for medical equipment
- + Delivered with power adapter suited for UK, EU, USA or AUS
- + When the battery pack is charged, the charger will continue as trickle charger

Charger 24V/DC SLA (MSP-C-AH-01100)

Specifications

Input voltage	100-240 VAC 50-60Hz
Input current	700-500mA
Output voltage	24 V/DC
Output current	400mA
IP class charger	IP40
Certification	UL2601-1 IEC60601-1 IEC60335-2-29
Operational temperature	0°C to +40°C
Storage temperature	-40°C to +70°C
Case colour	Black
Dimensions LxWxH	109x69x39mm

Order Information

Article number	MSP-C-AH-01101 for EU MSP-C-AH-01102 for GB MSP-C-AH-01103 for AUS MSP-C-AH-01104 for USA
Replacement for Arjo art. number	700-24201
Packaging	Individually packed, 50pcs per box
M.O.Q.	1 piece
Delivery Time	2 working days



Features

- + Suitable for Maxi Sky 600
- + Charger specifically designed for medical equipment
- + Delivered with power adapter suited for UK, EU, USA or AUS
- + When the battery pack is charged, the charger will continue as trickle charger